

ARTISAN
MEDICAL

Artisan Medical • 617 Stokes Road • Suite 4-306 • Medford, NJ 08055
ph: 866.930.3390 • fax: 609.939.0959

www.artisanmed.com



CABLE ORDERING INSTRUCTIONS

STEP 1	STEP 2	STEP 3	STEP 4	STEP 5	STEP 6	STEP 7
Determine Proximal Connector	Determine Distal Connector	Determine Fiber Optic Bundle	Determine Cable Length Size	Determine Instrument Connector Angle	Determine Sheathing Type • Non-Reinforced • 6" Reinforced • Fully-Reinforced	Determine Sheathing Color *(Fully-Reinforced offered in Gray only)*

STEP 1 PROXIMAL / ILLUMINATOR CONNECTOR

1  ACMI / CIRCON

2  WOLF / DYONICS

3  OLYMPUS / ACMI

4  STORZ

5  PILLING


STEP 2 DISTAL / INSTRUMENT CONNECTOR

6  WOLF SNAP ON

7  STORZ

8  ACMI / CIRCON SNAP ON

9  OLYMPUS

10  WOLF CLIP ON

11  PILLING

12  ACMI / CIRCON CLIP ON

STEP 3 FIBER OPTIC BUNDLE SIZE

Place gauge on bundle to determine diameter.



2.2 mm 3.0 mm 3.8 mm 4.8 mm 5.0 mm

STEP 4 FIBER OPTIC CABLE LENGTH

6 FT
8 FT
10 FT

Custom lengths available by request.

STEP 5 INSTRUMENT CONNECTOR ANGLE



Standard



45 Degrees

STEP 6 SHEATHING TYPE

0	**Non-Reinforced**
6	**6" Reinforced Ends**
FR	**Fully-Reinforced** *(Gray only)*

STEP 7 SHEATHING COLOR

CL	***Clear***	PK	***Pink***
WH	***White***	YL	***Yellow***
RD	***Red***	BK	***Black***
GY	***Gray***	GR	***Green***
BL	***Light Blue***		

Artisan Medical
617 Stokes Road • Suite 4-306
Medford, NJ 08055
ph: 866.930.3390
fax: 609.939.0959
www.artisanmed.com

MADE IN USA